AFIAS MxA/CRP

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
Figure 1. AFIAS MxA levels in viral infection(n=25), bacterial infection(n=39) and Non-infection(n=50) group.

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Clinical diagnosis	Sensitivity		Specificity	
	%	95% CI	%	95% CI
Viral infection	88.0	75.3 - 100.0	94.4	89.6 - 99.2
Bacterial infection	87.2	76.7 - 97.7	88.0	80.6 - 95.4

Table 1. Clinical sensitivity and specificity of AFIAS MxA/CRP to differentiate between viral infection and bacterial infection

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
- 1) Jonas, Olga B, et al, Drug-Resistant Infection: a threat to our economic future, World Bank Group, 2017
- 2) Clark D Russell et al., Co-infections, secondary infections, and antimicrobial use in patients hospitalized with COVID-19 during the first pandemic wave from the ISARIC WHO CCP-UK study: a multicentre, prospective cohort study, Lancet microbe, 2021
- 3) Antimicrobial Resistance Collaborators, Global burden of bacterial antimicrobial resistance in 2019: a systematic analysis, Lancet, 2022
- 4) AMR: Eroding the Foundations of Health Systems Across the EU and Around the World, Global Coalition on Aging,
- 5) Patrick Joseph and Eliot Godofsky, Outpatient Antibiotic Stewardship: A Growing Frontier—Combining Myxovirus Resistance Protein A With Other Biomarkers to Improve Antibiotic Use, OFID, 2022
- 6) Dan Coster et al., Using the kinetics of C-reactive protein response to improve the differential diagnosis between acute bacterial and viral infections, Infection, 2019

Boditech Med Inc.

43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gangwon-do, 24398, Korea

Tel: +82-33-243-1400 / Fax: +82-33-243-9373 / marketing@boditech.co.kr / www.boditech.co.kr/en